

SMARTdiagnostics 4.51 Web Release Notes

Expanded Media Support and a New Look for SMARTdiagnostics

This release introduces new personalization options and expanded media capabilities to improve how users interact with SMARTdiagnostics. Users can now explore an early beta version of Dark Mode, upload videos to support troubleshooting workflows, and interact with images more efficiently within the File Library.

New Features

- **Dark Mode Beta**

SMARTdiagnostics now includes a beta version of Dark Mode, introducing three new themes, Dark, Navy, and Night, in addition to the existing Light mode. Themes can be changed at any time from the color swatch menu in the upper-right corner of the application, giving users more flexibility to personalize their experience.

As this is an early beta feature, users may occasionally encounter minor visibility or display issues. If issues occur, users can switch back to Light mode at any time. Theme preferences will currently reset to Light mode after logging out.

- **Video Support in File Library**

Users can now upload videos to Issues, Assets, and the File Library, making it easier to capture and share machine behavior, troubleshooting context, and Motion Insights recordings from the mobile app.

Improvements to Existing Features

- **Enhanced Image Interactions in File Library**

Users can now click and drag images directly within the File Library for a smoother and more efficient image management experience when reviewing and organizing files.

Bug Fixes & General Improvements

This release includes additional fixes and performance enhancements to ensure a smooth and reliable experience.